


## Product Datasheet

### Pig IFN gamma ELISA Kit (orb315069)

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Pig IFN gamma ELISA Kit                                                                       |
| <b>Reactivity</b>          | Porcine                                                                                       |
| <b>Range</b>               | 15.6 pg/ml - 1000 pg/ml                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                         |
| <b>Target</b>              | IFNG                                                                                          |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                         |
| <b>Application notes</b>   | Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates                |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma, Tissue Homogenates                                    |
| <b>Sensitivity</b>         | 7 pg/ml                                                                                       |
| <b>Expiration Date</b>     | Please enquire.                                                                               |



Representative  
Standard  
Curve