


Product Datasheet

Human IGFBP4 ELISA Kit (orb315071)

Description	Human IGFBP4 ELISA Kit
Reactivity	Human
Range	500 pg/ml - 32,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D22-E258
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27934 MW
Uniprot ID	P22692
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 18.5-273 ng/ml.
Expiration Date	Please enquire.



Human
IGFBP4
PicoKine
ELISA Kit
standard...