


Product Datasheet

Mouse OSM/Oncostatin M ELISA Kit (orb315073)

Description	Mouse OSM/Oncostatin M ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-R206
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	30114 MW
Uniprot ID	P53347
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml.
Expiration Date	Please enquire.



Mouse
Oncostatin
M/OSM
PicoKine
ELISA Ki...