


Product Datasheet

Human SCF/KITLG/Kit ligand ELISA Kit (orb315074)

Description	Human SCF/KITLG/Kit ligand ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E26-A190
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	30899 MW
Uniprot ID	P21583
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 446-1729 pg/ml.
Expiration Date	Please enquire.



Human SCF
PicoKine
ELISA Kit
standard
CU...