


Product Datasheet

Human ADAM8/Ms2 ELISA Kit (orb315078)

Description	Human ADAM8/Ms2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-P497
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	88771 MW
Uniprot ID	P78325
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 33-220 pg/ml.
Expiration Date	Please enquire.



Human
ADAM8
PicoKine
ELISA Kit
standard ...