


Product Datasheet

Rat MIP-2 ELISA Kit (orb315081)

Description	Rat MIP-2 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S32-N100
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10783 MW
Uniprot ID	P30348
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Rat CXCL2
PicoKine
ELISA Kit
standard
cu...