


Product Datasheet

Mouse XCL1 / Lymphotactin / SCM 1 ELISA Kit (orb315088)

Description	Mouse XCL1 / Lymphotactin / SCM 1 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V22-G114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	12467 MW
Uniprot ID	P47993
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Mouse
XCL1/Lymphotactin
PicoKine ELISA K...