


Product Datasheet

Human KLK8/Kallikrein-8 ELISA Kit (orb315089)

Description	Human KLK8/Kallikrein-8 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q29-G260
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	28048 MW
Uniprot ID	O60259
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 600 pg/ml.
Expiration Date	Please enquire.



Human
KLK8/Kallikrein-
8 PicoKine
ELISA K...