


Product Datasheet

Rat CXCL5/LIX/ENA-78 ELISA Kit (orb315090)

Description	Rat CXCL5/LIX/ENA-78 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A38-Q130
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	14263 MW
Uniprot ID	P97885
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 20 ng/ml.
Expiration Date	Please enquire.



Rat CXCL5
PicoKine
ELISA Kit
standard
CU...