


Product Datasheet

Human HO-1/HMOX1/HSP32 ELISA Kit (orb315091)

Description	Human HO-1/HMOX1/HSP32 ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-M288
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	32819 MW
Uniprot ID	P09601
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 60-150 pg/ml.
Expiration Date	Please enquire.



Human HO-1/HMOX1 PicoKine ELISA Kit stan...