


Product Datasheet

Human AggreCAN/ACAN ELISA Kit (orb315093)

Description	Human AggreCAN/ACAN ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V20-R675
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and urine.
MW	250193 MW
Uniprot ID	P16112
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1.6 ng/ml.
Expiration Date	Please enquire.



Human
AggreCAN
PicoKine
ELISA Kit
standa...