


Product Datasheet

Human TNFRSF16 / NGFR / CD271 ELISA Kit (orb315095)

Description	Human TNFRSF16 / NGFR / CD271 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K29-N250
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	45183 MW
Uniprot ID	P08138
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml.
Expiration Date	Please enquire.



Human
TNFRSF16/NGFR
PicoKine ELISA
Kit s...