


Product Datasheet

Human S100A12/En Rage ELISA Kit (orb315101)

Description	Human S100A12/En Rage ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-E92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10575 MW
Uniprot ID	P80511
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:300, concentration in serum and plasma is around 80 ng/ml.
Expiration Date	Please enquire.



Human
S100A12
PicoKine
ELISA Kit
standar...