


Product Datasheet

Mouse beta IG-H3/TGFBI ELISA Kit (orb315104)

Description	Mouse beta IG-H3/TGFBI ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G24-H683
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	74597 MW
Uniprot ID	P82198
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is 200-626 ng/ml.
Expiration Date	Please enquire.



Mouse βIG-H3/TGFBI PicoKine ELISA Kit s...