


Product Datasheet

Mouse FGF23 ELISA Kit (orb315108)

Description	Mouse FGF23 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Y25-V251
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27758 MW
Uniprot ID	Q9EPC2
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 40-150 pg/ml.
Expiration Date	Please enquire.



Mouse
FGF23
PicoKine
ELISA Kit
standard ...