


Product Datasheet

Mouse FGF21 ELISA Kit (orb315110)

Description	Mouse FGF21 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A29-S210
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	23237 MW
Uniprot ID	Q9JJN1
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 500-700 pg/ml.
Expiration Date	Please enquire.



Mouse
FGF21
PicoKine
ELISA Kit
standard ...