


Product Datasheet

Human Serpin C1/Antithrombin-III ELISA Kit (orb315113)

Description	Human Serpin C1/Antithrombin-III ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H33-K464
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	52602 MW
Uniprot ID	P01008
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:20000, concentration in serum and plasma is 57-176 ug/ml.
Expiration Date	Please enquire.



Human Serpin
C1/Antithrombin-
III PicoKin...