


Product Datasheet

Mouse Epiregulin / EREG ELISA Kit (orb315115)

Description	Mouse Epiregulin / EREG ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V53-L101
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	18474 MW
Uniprot ID	Q61521
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 125 pg/ml.
Expiration Date	Please enquire.



Mouse
Epiregulin
PicoKine
ELISA Kit
stan...