


Product Datasheet

Human Complement Factor H/CFH ELISA Kit (orb315117)

Description	Human Complement Factor H/CFH ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from human plasma
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	139096 MW
Uniprot ID	P08603
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:20000, concentration in serum and plasma is around 300 ug/ml.
Expiration Date	Please enquire.



Human
Complement
H/CFH
PicoKine
ELISA Ki...