


## Product Datasheet

# Mouse Complement Factor H/CFH ELISA Kit (orb315118)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Complement Factor H/CFH ELISA Kit                                                                                                                                           |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                             |
| <b>Range</b>               | 0.78 ng/ml - 50 ng/ml                                                                                                                                                             |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: S875-V1252                                                                                                               |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).                                                                                                        |
| <b>MW</b>                  | 139138 MW                                                                                                                                                                         |
| <b>Uniprot ID</b>          | <b>P06909</b>                                                                                                                                                                     |
| <b>Sensitivity</b>         | 20 pg/ml                                                                                                                                                                          |
| <b>Dilution Range</b>      | Dilution ratio of 1:1000, concentration in serum and plasma is 5-20 ug/ml.                                                                                                        |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                   |



Mouse  
Complement  
H/CFH  
PicoKine  
ELISA Ki...