


Product Datasheet

Mouse Cripto/TDGF1 ELISA Kit (orb315122)

Description	Mouse Cripto/TDGF1 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: R26-Q150
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	18754 MW
Uniprot ID	P51865
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 62.5 pg/ml.
Expiration Date	Please enquire.



Mouse
Cripto/TDGF1
PicoKine
ELISA Kit st...