

Product Datasheet

IL8/CXCL8 Antibody (orb315129)

Description

IL8/CXCL8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E. coli-derived human IL-8 recombinant protein (Position: A23-S99).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

11098 MW

Uniprot ID


P10145

Dilution Range


Western blot, 0.1-0.5µg/ml, Human
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat ELISA , 0.1-0.5µg/ml, Human, -

Expiration Date

12 months from date of receipt.



WB analysis of IL-8 using anti-IL-8 antibody.



IHC analysis of IL-8 using anti-IL-8 antibody.