

Product Datasheet

Gremlin 1/GREM1 Antibody (orb315134)

Description

Gremlin 1/GREM1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human Gremlin 1 (151-184aa TMMVTLNCPQLQPPTKKRVTRVKQCRCISIDLD), identical to the related mouse and rat sequences.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

20697 MW

Uniprot ID

[O60565](#)

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat

Expiration Date

12 months from date of receipt.

130KD —

70KD —

55KD —

35KD —

25KD —

15KD —

WB analysis
of Gremlin 1
using anti-
Grem...