



Product Datasheet

ANKS1B antibody (orb31515)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ANKS1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	138066
Uniprot ID	Q7Z6G8
NCBI	NP_001190994.1 , NP_690001.3 , NP_001190996.1 , NP_858056.2 , NP_001191009.1 , NP_001190995.1 , NP_001191010.1 , NP_001190999.1 , NP_001190998.1 , NP_001190997.1 , NP_064525.1 , NP_001191008.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of NCI-H292 cell lysate.