

## Product Datasheet

### Cyp3a13 antibody (FITC) (orb31519)

|                        |                                                                  |
|------------------------|------------------------------------------------------------------|
| <b>Description</b>     | Cyp3a13 antibody (FITC)                                          |
| <b>Species/Host</b>    | Rabbit                                                           |
| <b>Reactivity</b>      | Mouse                                                            |
| <b>Conjugation</b>     | FITC                                                             |
| <b>Immunogen</b>       | Recombinant Mouse Cytochrome P450 3A13 protein (1-503AA)         |
| <b>Target</b>          | Cyp3a13                                                          |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b> | Liquid                                                           |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>            | For research use only                                            |
| <b>Isotype</b>         | IgG                                                              |
| <b>Clonality</b>       | Polyclonal                                                       |
| <b>Uniprot ID</b>      | <a href="#">Q64464</a>                                           |
| <b>Dilution Range</b>  | ELISA                                                            |
| <b>Expiration Date</b> | 12 months from date of receipt.                                  |