

Product Datasheet

IL5 antibody (orb315196)

Description	Mouse monoclonal antibody to IL5
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, Immunostaining, WB
Target	IL-5
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Concentration	0.5 mg/vial
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Application notes	In vitro using serum free medium
Isotype	IgG1
Clonality	Monoclonal
NCBI	NP_000870
Entrez	3567
Expiration Date	12 months from date of receipt.