

Product Datasheet

IL17A antibody (orb315205)

Description	Mouse monoclonal antibody to IL17A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	IL-17A
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Entrez	3605
Expiration Date	12 months from date of receipt.