

Product Datasheet

IP-10 antibody (orb315213)

Description	Mouse monoclonal antibody to IP-10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	IP-10
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
NCBI	NP_001556
Entrez	3627
Expiration Date	12 months from date of receipt.