

Product Datasheet

MCP-1 antibody (orb315217)

Description	Mouse monoclonal antibody to MCP-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	MCP-1
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
NCBI	NP_002973
Entrez	6347
Expiration Date	12 months from date of receipt.