

Product Datasheet

IL10 antibody (orb315232)

Description	Mouse monoclonal antibody to IL10
Species/Host	Mouse
Reactivity	Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Target	IL-10
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Concentration	0.5 mg/vial
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Application notes	In vitro using serum free medium.
Isotype	IgG1
Clonality	Monoclonal
NCBI	NP_001038192
Entrez	694931
Expiration Date	12 months from date of receipt.