

## Product Datasheet

### IFN-gamma antibody (orb315237)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to IFN-gamma                               |
| <b>Species/Host</b>        | Mouse                                                                |
| <b>Reactivity</b>          | Rat                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                         |
| <b>Tested Applications</b> | ELISA, IHC-P, WB                                                     |
| <b>Target</b>              | IFN-gamma                                                            |
| <b>Preservatives</b>       | Product (vacuum-packed); Reconstitution with sterile distilled water |
| <b>Form/Appearance</b>     | Lyophilized                                                          |
| <b>Concentration</b>       | 0.5 mg/vial                                                          |
| <b>Storage</b>             | Store at 4°C until expiry date (mentioned on vial)                   |
| <b>Note</b>                | For research use only                                                |
| <b>Application notes</b>   | In vitro using serum free medium.                                    |
| <b>Isotype</b>             | IgG1                                                                 |
| <b>Clonality</b>           | Monoclonal                                                           |
| <b>NCBI</b>                | <b>NP_620235</b>                                                     |
| <b>Entrez</b>              | <b>25712</b>                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                      |