



## Product Datasheet

### MAGED4B antibody (orb31525)

|                            |                                                                                                                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Peptide affinity purified rabbit polyclonal antibody.                                                                                                                                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                 |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                                                                                                            |
| <b>Target</b>              | MAGED4                                                                                                                                                                                                                                                       |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.                                                                                 |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                   |
| <b>MW</b>                  | 81378                                                                                                                                                                                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q96JG8</a>                                                                                                                                                                                                                                       |
| <b>NCBI</b>                | <a href="#">NP_803881.1</a> , <a href="#">NP_001258990.1</a> , <a href="#">NP_001258991.1</a> , <a href="#">NP_001258992.1</a> , <a href="#">NP_001229291.1</a> , <a href="#">NP_110428.2</a> , <a href="#">NP_803879.1</a> , <a href="#">NP_001092270.1</a> |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:10-50                                                                                                                                                                                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                              |



Western blot analysis of HepG2 cell line...