

Product Datasheet

IL2 antibody (orb315250)

Description	Rat monoclonal antibody to IL2
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP, NeA, WB
Target	IL-2
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Concentration	0.5 mg/vial
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Application notes	In vitro using Hybridoma-SFM medium + 1% FCS
Isotype	IgG2b
Clonality	Monoclonal
NCBI	NP_032392
Entrez	16183
Expiration Date	12 months from date of receipt.