

Product Datasheet

IL4 antibody (orb315251)

Description	Rat monoclonal antibody to IL4
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL-4
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Concentration	0.5 mg/vial
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Application notes	In vitro using serum free medium
Isotype	IgG1
Clonality	Monoclonal
NCBI	NP_067258
Entrez	16189
Expiration Date	12 months from date of receipt.