

Product Datasheet

IL10 antibody (orb315255)

Description	Rat monoclonal antibody to IL10
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Target	IL-10
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
NCBI	NP_034678
Entrez	16153
Expiration Date	12 months from date of receipt.