

Product Datasheet

IFN-gamma antibody (FITC) (orb315260)

Description	Mouse monoclonal antibody against IFN-gamma conjugated to FITC
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	FITC
Tested Applications	FC, IHC
Target	IFN-gamma
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Concentration	50 µg/vial
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Application notes	In vitro using serum free medium
Isotype	IgG1
Clonality	Monoclonal
NCBI	NP_000610
Entrez	3458
Expiration Date	12 months from date of receipt.