

Product Datasheet

CD3 antibody (FITC) (orb315265)

Description	Mouse monoclonal antibody against CD3 conjugated to FITC
Species/Host	Mouse
Reactivity	Monkey
Conjugation	FITC
Tested Applications	FACS, WB
Target	CD3
Preservatives	Product (vacuum-packed); Reconstitution with sterile distilled water
Form/Appearance	Lyophilized
Concentration	50 µg/vial
Storage	Store at 4°C until expiry date (mentioned on vial)
Note	For research use only
Application notes	In vitro using serum free medium.
Isotype	IgG1
Clonality	Monoclonal
Entrez	705270, 699582, 699467
Expiration Date	12 months from date of receipt.