


Product Datasheet

G007-LK (orb315315)

Description	G007-LK is a potent, "rule of 5" compliant and a metabolically stable TNKS1/2 inhibitor.
Note	For research use only
Formula	C ₂₅ H ₁₆ CLN ₇ O ₃ S
MW	529.9576
CAS Number	[1380672-07-0]
SMILES	<chem>C1C=C([H])C([H])=C([H])C([H])=C1N1C(/C(/[H])=C([H])/C2=NN=C(C3C([H])=C([H])C(C#N)=C([H])C=3[H])O2)=NN=C1C1C([H])=C([H])C(=C([H])N=1)S(C([H])([H])([H])(=O)=O</chem>
Expiration Date	12 months from date of receipt.



Chemical
structure of
G007-LK