

Product Datasheet

PI3K-IN-2 (orb315427)

Description	PI3K-IN-2 is a PI3K inhibitor, inhibits pPKB and pS6 with "++++/+++(+)" in A2058 melanoma cell.
Note	For research use only
Formula	$C_{17}H_{20}F_3N_7O_2$
MW	411.38
CAS Number	[1225037-39-7]
SMILES	<chem>C(N)1=NC=C(C2=NC(N3CCOCC3)=NC(N3CCOCC3)=N2)C(C(F)(F)F)=C1</chem>
Expiration Date	12 months from date of receipt.