

Product Datasheet

MID1 antibody (orb31543)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MID1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75251
Uniprot ID	O15344
NCBI	NP_150632.1 , NP_000372.1 , NP_001180209.1 , NP_001180210.1 , NP_001180206.1 , NP_001180208.1 , NP_001092094.1 , NP_001180207.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

