


Product Datasheet

ORM1 protein (orb315543)

Description	ORM1 protein
Species/Host	E. coli
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	ORM1
Form/Appearance	20mM Tris-HCl, 50mM NaCl, 2M Urea, 1mM DTT, pH 8.0
Concentration	0.4 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Application notes	N-His6
Protein Sequence	QIPLCANLVP VPITNATLDR ITGKWFYIAS AFRNEEYNKS VQEIQATFFY FTPNKTEDTI FLREYQTRQD QCIYNTTYLN VQRENGTISR YVGGQEHFAH LLILRDTKTY MLAFDVNDEK NWGLSVYADK PETTKEQLGE FYEALDCLRI PKSDVVYTDW KKDKCEPLEK QHEKERKQEE GES
Purity	> 90.0% as determined by 12% SDS-PAGE
MW	23.9 kDa
NCBI	000607
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
ORM1 protein