



Product Datasheet

IL1RA antibody (orb315548)


Description	IL1RA antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human IL1RA . Please contact us for the exact immunogen sequence. The peptide is available as orb374902.
Target	IL1RA
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 100 µl): 1 mg/ml- 200 µg (in 200 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P25085, P25086, P18510
NCBI	000577, 000568
Entrez	3557
Dilution Range	IHC-P: 1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat panc...



IHC-P staining of rat colon tissue using...