

Product Datasheet

NKp44 antibody (orb315559)

Description

NKp44 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human NKp44. Please contact us for the exact immunogen sequence. The peptide is available as orb374913.

Target

NKp44

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 100 µl): 1 mg/ml- 200 µg (in 200 µl): 1 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

Source

The anti-GFP antibody is a rabbit polyclonal antiserum obtained from rabbits immunised with highly purified native GFP derived from Aequorea victoria jellyfish.

Uniprot ID

[O95944](#)

NCBI

[21](#), [31](#), [004819](#), [004828](#)

Entrez


[9436](#)

Dilution Range


WB: 1:100-800, IHC-P: 1:100-500

Expiration Date


12 months from date of receipt.



Western blot analysis of mouse lymph node...



Immunohistochemical staining of human endometrial tissue...



IHC-P image of mouse lymph node tissue u...