



Product Datasheet

Kir3.1 (Phospho-S185) antibody (orb315577)


Description	Rabbit polyclonal antibody to Kir3.1 (Phospho-S185)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S185 of human Kir3.1 protein. The exact sequence is proprietary.
Target	KCNJ3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P48549, P63251, P63250
Entrez	50599, 16519, 3760
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Kir3.1 (Phospho...



Immunohistochemical analysis of Kir3.1 (...)



Immunofluorescent analysis of Kir3.1 (Ph...