


Product Datasheet

Cytokeratin 8 (Phospho-S432) antibody (orb315586)

Description

Rabbit polyclonal antibody to Cytokeratin 8



Western blot analysis of Cytokeratin 8 (...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S432 of human Cytokeratin 8 protein. The exact sequence is proprietary.

Target

KRT8

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source


Rabbit

Uniprot ID
[Q10758](#), [P11679](#), [P05787](#)
Entrez
[3856](#), [16691](#), [25626](#)
Dilution Range


WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500

Expiration Date

12 months from date of receipt.



Immunohistochemical analysis of Cytokera...



Immunofluorescent analysis of Cytokera...