



Product Datasheet

LCK (Phospho-Y505) antibody (orb315588)

Description	Rabbit polyclonal antibody to LCK (Phospho-Y505)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y505 of human LCK protein. The exact sequence is proprietary.
Target	LCK
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P06240, Q01621, P06239
Entrez	313050, 3932, 16818
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200
Expiration Date	12 months from date of receipt.



Western blot analysis of LCK (Phospho-Y505) antibody.



Immunohistochemical analysis of LCK (Phospho-Y505) antibody.