



Product Datasheet

LYN (Phospho-Y397) antibody (orb315590)

Description	Rabbit polyclonal antibody to LYN (Phospho-Y397) also known as p53Lyn and p56Lyn whose expression ca...
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y397 of human LYN protein. The exact sequence is proprietary.
Target	LYN
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P25911, Q07014, P07948
Entrez	17096, 4067, 81515
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500
Expiration Date	12 months from date of receipt.



Western blot analysis of LYN (Phospho-Y3...



Immunohistochemical analysis of LYN (Pho...



Immunofluorescent analysis of LYN (Phosph...