

Product Datasheet

JNK1/2/3 (Phospho-T183/Y185) antibody (orb315601)

Description

Rabbit polyclonal antibody to JNK1/2/3

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IF, WB

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T183/Y185 of human JNK1/2/3 protein. The exact sequence is proprietary.

Target

MAPK8; MAPK9; MAPK10

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[Q9WTU6](#), [P49186](#), [P49185](#), [P45984](#), [P49187](#), [P53779](#), [Q91Y86](#), [P45983](#), [Q61831](#)

Entrez


[116554](#), [25272](#), [26419](#), [26420](#), [5602](#), [5599](#), [50658](#), [5601](#)

Dilution Range


WB: 1-500:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of JNK1/2/3 (Phospho...)



Immunofluorescent analysis of JNK1/2/3 (...)