



Product Datasheet

MEF2A (Phospho-T312) antibody (orb315611)

Description	Rabbit polyclonal antibody to MEF2A (Phospho-T312)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T312 of human MEF2A protein. The exact sequence is proprietary.
Target	MEF2A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q2MJT0, Q60929, Q02078
Entrez	309957, 4205, 17258
Dilution Range	WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500, IP: 1-10:100
Expiration Date	12 months from date of receipt.



Western blot analysis of MEF2A (Phospho-...



Immunohistochemical analysis of MEF2A (P...