


Product Datasheet


Ezrin/Radixin/Moesin (Phospho-T567/564/558) antibody (orb315620)

Description


Rabbit polyclonal antibody to



Western blot analysis
of Ezrin/Radixin/M...



Immunohistochemical
analysis of Ezrin/Ra...



Immunofluorescent
analysis of
Ezrin/Radi...

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Primate, Rat,
Zebrafish

Conjugation

Unconjugated

Tested Applications

IF, IH, WB

Immunogen

KLH-conjugated synthetic phosphopeptide
corresponding to residues surrounding
T567/564/558 of human
Ezrin/Radixin/Moesin protein. The exact
sequence is proprietary.

Target

MSN; RDX; EZR

Preservatives

Liquid in 0.42% Potassium phosphate,
0.87% Sodium chloride, pH 7.3, 30%
glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate,
0.87% Sodium chloride, pH 7.3, 30%
glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and
store at -20°C for one year. Avoid
freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

P31977, **P35241**, **P15311**, **P26038**,
P26041, **O35763**, **P26040**, **P26043**

Entrez

4478, **17698**, **5962**, **19684**, **54319**,
22350, **81521**

Dilution Range

WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC:
1-100:500

Expiration Date

12 months from date of receipt.